


Dental Gelatin and Collagen Cubes

Absorbable Micro-Porous Hemostatic Cubes

Dental Bleeding Is Now Under Control



Ordering Information

Code	Description	Device Specifications	Quantity (per box)
201CD001	Colident™ Dental Cubes	10 x 10 x 10mm	60 pieces
201CD002	Colident™ Dental Cubes	10 x 10 x 10mm	40 pieces
201CD003	Colident™ Dental Cubes	10 x 10 x 10mm	20 pieces

incura

To know more about our products,
contact us via:

website

www.incura-med.com

mobile

+966 54 9378601
+966 54 9384285
+201 229 050709
+201 028 643738

email

wesam.sarhan@incura-med.com
mousa.elsayed@incura-med.com
info@incura-med.com

address

7775 St. 241, AlYasmin District,
Riyadh 13329, Kingdom of
Saudi Arabia



incura




Advanced Hemostatic Solution for Dental Applications


Colident™ enhances blood clotting and platelet aggregation with an absorption capacity of 55 times its original weight to provide dentists with an efficient hemostasis tool they can depend on.

Colident™'s biodegradable structure allows for:

- Instant blood absorption and bleeding control
- High blood absorption capacity
- Accelerated hemostasis and platelet aggregation




Mechanism of Action




1 Immediate Absorption

Upon blood contact, the highly porous structure and the hydrophilic nature of Colident™ sponge allow for fluid absorption, without the need for pre-wetting




2 Platelet Aggregation

As blood is absorbed, platelets and clotting factors are concentrated within the sponge bleeding site interface.



3 Clot Formation

Concentrated clotting factors accelerate the clotting process, helping to seal the bleeding site and reduce blood loss while preserving structural stability.



4 Biodegradation

Colident™ is fully biodegradable, breaking down into harmless byproducts that are absorbed by the body in less than 10 days, leaving no residue at the surgical site.


Colident™ 70X magnification

Immediate Action

Colident™ sponges feature a hydrophilic, cross-linked porous structure, allowing for fast, efficient absorption without prewetting during critical surgical procedures

Super Absorption Capacity

Colident™ sponges are meticulously engineered with a super-absorbent, porous structure, absorbing blood three times more effectively than Ultrafoam, Cutanplast, and Curaspon, ensuring superior bleeding control and superior mechanical stability



Fully Biodegradable

Colident™ sponges, made from halal gelatin, are fully biodegradable and absorbed by the body in less than 10 days, supporting natural healing